

The Intertek logo is located in the top right corner. It features the word "intertek" in a lowercase, sans-serif font. The "i" is a small dot, and the "t"s have a unique, slightly curved top. The color is a bright yellow.

Total Quality. Assured.

INTERTEK PHARMACEUTICAL SERVICES

BIPHARMACEUTICAL SERVICES

Expert Analytical & Bioanalytical Services



Our Biopharmaceutical Expertise

To meet your development milestones and comply with regulatory requirements, you will need high quality bioanalytical and analytical data to aid informed decision-making and identify sources of risk. Safety assessment is key, as slight changes in the structure, physicochemical properties, stability and impurity profile of a biologic can provoke an adverse immune response.

Our thought-leaders have nearly 30 years of experience in biopharmaceutical development support across a wide range of product types. We provide regulatory-led, phase-appropriate, tailored analytical program design with GLP / GCP / GMP compliant laboratory services which help you to navigate the challenges of development, regulatory submission, and manufacturing.

With facilities located in India, Europe (UK, France, Switzerland) and the USA, our strong scientific and technical leadership coupled with project management and regulatory support, deliver a responsive, flexible (bio)analytical resource, to drive your development and manufacturing programs forward.

Our biopharmaceutical services

- GLP Clinical & Preclinical Bioanalysis (PK, ADA, Nab)
- Immunogenicity Studies
- Analytical Programme Design in line with ICH Q6B
- Inhaled & Nasal Drug Development
- Structural Characterisation
- Physicochemical Properties
- Biophysical Characterisation
- Comparability Studies
- Biosimilar Programmes
- Process Residuals Determination
- Product Related Impurities Determination
- Purity and Impurity Assessment
- GMP Potency / Cell Based Assays
- Method Development & Validation
- Extractables / Leachables
- GMP Quality Control Testing
- GMP Batch & Final Product Release Testing
- ICH Stability Studies
- Forced Degradation Studies

Our experience

- Recombinant Proteins & Glycoproteins
- Monoclonal Antibodies
- Antibody-Drug Conjugates
- Biosimilars / Biobetters
- Peptides
- Growth Factors
- Interferons
- Interleukins
- Vaccines
- Viral Vectors
- Oligonucleotide Therapeutics

To achieve a well-characterised biologic, we apply a strategic analytical program that is tailored to your molecule

Protein Characterisation

Protein therapeutics can present significant analytical challenges. Our characterisation packages are tailored to your biomolecule and reflect the requirements of ICH Guidelines to meet the analytical challenges of a 'well characterized' or 'specified' biological product.

Programs encompass many different analytical techniques and provide information ranging from evaluation of physicochemical properties and structural features including primary, secondary and higher order structure and assessment of post-translational modifications. Programs also include determination of biological potency and assessment of purity /impurity profiles. We continually invest in advanced analytical instrumentation which allows us to deliver data with the highest sensitivity, accuracy and resolution.



Structural Characterisation	
Amino Acid Sequencing/ Peptide Mapping:	Sequencing studies and peptide mapping using a broad range of enzymatic or chemical digestion followed by Mass Spectrometry (LC-MS/MS or MALDI-TOF MS).
Amino acid composition:	Pharmacopeia methods.
Terminal amino acid sequence:	Confirmation of N- and C-terminal sequences and evaluation of modifications and / or heterogeneity.
Disulphide bridge mapping:	Assessment of the degree and positions of both expected and mis-matched disulphide bridges by extended LC-MS/MS study and colorimetric test for free sulphhydryl groups.
Carbohydrate structure:	Glycosylation studies typically including levels of monosaccharides and sialic acid, N/O linked glycan profiling (NPLC, HILIC, IEX or CE-LIF), enzymatic digest and MALDI -TOF or LC-MS/MS.
Physico-chemical Properties	
Molecular weight:	Molecular weight of intact proteins by MALDI-MS, ESI-MS and LC-MS supported by orthogonal techniques such as MALLS and SDS-PAGE.
Isoform pattern:	Isoform and impurity studies using PAGE, SDS-PAGE, IEF, CE, HPLC.
Extinction coefficient:	Determination and Validated Extinction Coefficient studies.
LC patterns:	For ID, homogeneity, purity – HPLC, UPLC, SEC, RP HPLC, IEX, AEX.
Spectroscopic patterns:	CD, FTIR, 1D & 2D NMR, Fluorescence, UV-Visible.
Electrophoretic patterns:	CE(CZE), cIEF, CGE, SDS and NATIVE PAGE, Western Blot.
Concentration:	Lowry, BCA, Total AA, Total Nitrogen, Bradford.
Aggregation studies:	SEC (MALLS), DLS, AUC, Western Blot, CE, Gel Electrophoresis, SEM/TEM. Cryo-TEM.
Process Impurities:	Residual host cell DNA by qPCR, Residuals (such as antibiotics, antifoaming agents).
Potency Assays	
Cell-based Assays:	To support characterisation, stability, comparability testing and product release, for example, Complement-dependent cytotoxicity (CDC), Antibody-dependent cell cytotoxicity (ADCC) and Neutralisation and Proliferation Assays.

We apply our bioanalytical expertise and industry insight to design strategic and efficient bioanalytical programs

Bioanalytical Expertise

Large molecule bioanalysis

We have extensive experience in the development, validation, and sample analysis of quantitative and qualitative Good Laboratory Practice (GLP) and Good Clinical Practice (GCP) immunochemistry and LC-MS/MS expertise in support of pre-clinical and clinical development programs in full regulatory compliance (FDA, EMA and OECD GLP regulatory standards).

Our bioanalytical experts are highly experienced in developing and validating assays for pharmacokinetic (PK), toxicokinetic (TK), pharmacodynamic (PD), immunogenicity, efficacy and safety biomarker studies with diverse immunochemistry technologies, radioimmunochemistry (RIA) and functional assay platforms.

We offer comprehensive bioanalytical services to support complex products including proteins, synthetic peptides, humanized monoclonal antibodies, antibody drug conjugates, conjugated drugs, growth factors, hormones and cytokines. Additionally, we offer industry unique capabilities, for example, quantitative, high resolution NMR which is an ideal tool for bioanalysis of PEGylated moieties.

- Quantitative Ligand Binding Assay Capabilities
- Quantitative ELISAs for Proprietary Compounds
- Immunogenicity Studies
- Neutralization Cell-based Assay Development, Validation and Sample Analysis
- Enzymatic Assays
- Fluorometric Assays, Luminescence Assays
- Mode of Action Studies
- Bioanalytical LC-MS/MS for Biologics
- NMR approaches for PEGylated Biomolecules
- Biomarker Assays

Biomarker assays

Intertek's GLP/GCP/GMP compliant laboratories provides support to clients focused on the development of pharmaceuticals and biological medicines. Our biomarker services provide fit-for-purpose validations in support of your context of use. We provide a continuous and iterative process which evolves with your context of use for the biomarker data.

Intertek deploy a diverse array of analytical platforms to support biomarker analysis including ELISA, flow cytometry, cell-based assays, ECL multiplex platforms and high sensitivity measurements by the Quanterix Simoa HD-X Analyzer™.

CASE STUDY

Novel Approaches for An Enzyme Activity Assay

A client desired an activity assay for a PEGylated enzyme for which a commercial colorimetric assay was available. The colorimetric assay did not meet the performance criteria for regulated work. An immunoassay was then developed but was subject to significant matrix effects.

Our Solution

Development of a replacement assay was complicated by endogenous substrate and enzymatic product. To overcome this, the specificity of LC-MS/MS was employed and an activity assay developed using a stable labeled substrate which produced a labeled product, which could be differentiated from the endogenous analyte.

Benefit Delivered to our Client

An enzymatic activity LC-MS/MS assay was developed and successfully validated to regulatory standards. The method was used in multiple pharmacokinetic studies and enabled the client to move forward with their drug development program.

Specialist Expertise

Biosimilar services

Our biosimilar thought leaders apply a strategic approach to developing biosimilarity data with programs that provide highly relevant early stage characterisation and late stage comparative data. These programs begin with extensive structural and functional characterization of both the proposed biosimilar product and the reference product and evaluate and compare all pertinent features of the biosimilar product and are based on the criteria outlined in ICH Q6B. This analytical characterization serves as the foundation of a biosimilar development program. Programs encompass many different analytical techniques and provide information ranging from evaluation of physicochemical properties and structural features including primary, secondary and higher order structure and assessment of post-translational modifications and cell-based assays to determination of biological potency with the goal to show that the molecule is structurally and functionally highly similar to a reference product and so is anticipated to behave like the reference product.

Monoclonal antibodies

Our experience spans recombinant monoclonal antibodies and related products such as biosimilars, fusion proteins, Fab-fragments and Fc fragments and antibody drug conjugates (ADCs) throughout the product lifecycle with a focus on monitoring relevant critical quality attributes (CQAs), assessment of the impact of process changes on physicochemical properties and structure, the presence of product-related impurities or process-related impurities.

- Characterisation (Q6B)
- Primary Sequence Confirmation (Peptide Map)
- N/C terminal Sequence, Disulphide Bridge Mapping and PTMs
- Carbohydrate Structure
- Aggregation Studies
- Higher Order Structural Studies
- Aggregation Studies
- Heterogeneity, Purity and Variants
- Potency

Protein Analytics

Post Translational Modifications (PTM)

We apply our knowledge of anticipated PTM and potential PTMs based on your product's specific structural characteristics to design strategic analytical programs that meet the requirements of regulatory expectations and industry guidance such as ICH Q6B. Using enzyme or chemical digestion coupled with highly sensitive liquid chromatography tandem mass spectrometry (LC-MS/MS) and a range of chromatographic techniques we can detect and identify a range of PTMs.

- Deamidation
- Glycosylation
- Phosphorylation
- Acetylation
- Alkylation
- Sulfation
- Glycylation
- Methylation
- Oxidation
- Mis-matched S-S bridges
- Truncation
- N/C-terminal Modifications

Biophysical characterisation

Our biophysical characterisation team applies a wide range of techniques such as to interrogate the biophysical behavior of your molecule. We examine product's higher-order structure (HOS), including secondary and tertiary structure using Circular Dichroism (CD), FTIR, Fluorescence, 1D / 2D High field 600MHz NMR studies (NOESY & TOCSY), aggregation (using SEC-MALS, DLS, SV-AUC), oligomerization and degradation of the drug substance. To support formulation development we study proteins in the buffer and in the presence of additives / excipients of interest. To support IND or NDA submissions we assess spectroscopic, thermodynamic and hydrodynamic properties.

Potency assays

Potency bioassays are an integral part of release and stability testing. Intertek's experienced specialists perform cell based potency and ligand binding assays which can be performed to GMP for regulatory submissions or as a non-regulatory study under high quality standards.

Glycosylation

Glycosylation studies are designed to be product specific, however, these typically include determination of the levels of neutral and amino monosaccharides as well as sialic acids, assessment of glycoform distribution and glycan structure elucidation. Multiple technologies are applied to the determination including selective enzymatic cleavage and MALDI-TOF Mass spectrometry HPLC, HILIC, IEX or CE-LIF, to provide the level of structural information required.

Product related impurities determination


Our scientists perform detailed characterization using a diverse range of technologies which include chromatographic, electrophoretic, mass spectrometry (MALDI-MS, LC-MS/MS) and spectroscopic (IR, NMR and fluorescence approaches) to determine truncated forms and other modified forms such as deamidated, isomerised, mismatched S-S linked, oxidised or altered conjugated forms (e.g. glycosylation, phosphorylation). For truncated forms, we use an enzyme or chemical cleavage of peptide followed by characterisation by chromatography (e.g. HPLC) or SDS-PAGE approaches. We also perform peptide mapping to yield useful information about the variant.

Process residuals determination

Our process residual analysis services apply a range of highly specific and sensitive techniques to residual cell substrates (e.g. host cell proteins, host cell DNA), cell culture (e.g. inducers, antibiotics, or media components), or chromatographic media used in purification, solvents and buffer components to offer robust quantification supporting process validation, monitor the batch to batch variation and support GMP lot release. For trace antibiotics such as Kanamycin, Tobramycin, Gentamycin, Amoxicillin, Chloramphenicol, we apply highly selective approaches such as tandem mass spectrometry (LC-MS/MS) with Multiple Reaction Monitoring (MRM). We also apply MS techniques to the determination of surfactants and antifoams agents which typically include Triton-X, Tween 20, Tween 80, Pluronics and polyglycol P2000.

- Kanamycin
- Tobramycin
- Gentamycin
- Amoxicillin
- Chloramphenicol
- Methotrexate
- Glutathione
- Dithiothreitol
- IPTG
- Triton-X
- Tween
- Pluronics
- Polyglycol P2000
- Trace metals
- Downstream-derived impurities
- Leachables from process equipment
- Particulates

Development Support



ICH stability testing

With dedicated ICH storage stability facilities that are integrated with advanced analytical capabilities, we can offer stability-indicating method development and validation alongside efficient stability programs that focus on the critical quality attributes of your molecule.

Inhaled & intranasal biologic drug development

Formulation of biopharmaceuticals for inhaled or nasal delivery allows for a more convenient method of administering compounds systemically, and also allows direct targeting of the respiratory system, however, this brings new challenges. Our experts integrate formulation with tailored analytics that focuses on the impact of formulation on the safety and efficacy of the drug whilst monitoring aggregation, fragmentation, deamidation, hydrolysis, oxidation, isomerisation, succinimidation,

deglycosylation, disulphide bond formation/breakage and other crosslinking reactions.

Our Intertek Centre of Excellence for Inhaled Biologics deploys a strategic programme of formulation, analytical testing, and clinical manufacturing services. Our orthogonal analytical approach aims to fully characterise the biologic drug substance whilst assessing whether the device delivery mechanism has adversely affected the product performance.

- CMC Development & Analytical Support
- Biologics Characterisation (Q6B)
- Formulation Development
- Method Development/Validation
- Device Screening
- ICH & Accelerated Stability Studies
- Product Characterisation Studies



Americas

 USA +1 800 967 5352

Europe

UK +44 161 721 5247

Asia

India +91 9819117214
+91 7045312819
China +86 21 6073 7735

 biopharma@intertek.com

 intertek.com/pharmaceutical